



Practitioner

Rosemary Hill, BSN, CWOON, CETN(C), with over 12 years of wound care experience, oversees wound care for inpatients and outpatients at Lions Gate Hospital, Vancouver Coastal Health, located in North Vancouver, Canada. Rosemary is currently President of the Canadian Association for Enterostomal Therapy (CAET).



Patient Condition

74-year-old female with a venous leg ulcer for longer than 5 months and confirmed hemangioma. The patient had previously received multiple courses of oral antibiotics, with minimal effect on wound progression, prior to her referral to the wound specialty team at Lions Gate Hospital, Vancouver Coastal Health for stalled wound healing.

Point of Care Identification of *Pseudomonas aeruginosa* with the MolecuLight *i:X*[®] Led to Targeted Treatment

Pseudomonas aeruginosa is one of the most common pathogens found in chronic leg ulcers, frequently leading to a stall in wound healing¹⁻². Its unique intrinsic and acquired antibiotic resistances make early identification and the selection of *Pseudomonas*-specific treatment regimes critical in wound care¹⁻². *P. aeruginosa* sometimes presents with unique clinical symptoms, including a malodorous, greenish crust and a greenish tinge on removed dressings³, but there are cases when no specific signs and symptoms are observed³. This bacterial species has unique fluorescent properties and is indicated by a cyan color on MolecuLight *i:X* fluorescence images⁴, enabling immediate species identification.

This patient's stalled wound exhibited no clinical signs and symptoms specific to *Pseudomonas aeruginosa* colonization, yet real-time visualization of cyan fluorescence on MolecuLight *i:X* images (Figure 2) strongly suggested *P. aeruginosa* (>10⁴ CFU/g)⁴. Obtaining this information at the point of care led the clinician to immediately select an antimicrobial dressing indicated for use against *P. aeruginosa*. Swabs taken from regions of cyan fluorescence under MolecuLight *i:X* fluorescence guidance later confirmed moderate growth of *P. aeruginosa*.



Figure 1: Standard Image. Image taken under standard lighting conditions.




Figure 2: Fluorescence Image. Arrows indicate regions of cyan fluorescence. Swabs of cyan regions later confirmed to be moderate growth of *P. aeruginosa*.

At the next visit, cyan fluorescence was no longer detected, supporting the effectiveness of the chosen antimicrobial treatment. Systemic antibiotics were not re-prescribed. With this effective eradication of *P. aeruginosa* contamination, the wound began to heal within three weeks.

MolecuLight *i:X*[®]

The MolecuLight *i:X* allows clinicians to quickly, safely and easily identify wounds with bacteria⁴⁻⁷ (at loads of $>10^4$ CFU/g, in combination with CSS) and measure wounds^{5,7} at the point of care to provide them with valuable information to inform treatment and monitor progress^{6,7}.

Testimonial

"I would not have suspected *Pseudomonas* on this wound. Both the Infectious Disease physician and myself were surprised when the *i:X* images showed clear presence of the cyan color, which influenced my clinical practice."



— Rosemary Hill, BSN, CWOCN, CETN(C)

Visit www.moleculight.com

+1.647.362.4684

Toll Free 1.877.818.4360 (Canada)

info@moleculight.com

Follow us:  



References:

1. Serra R *et al.* Expert Rev Anti Infect 2015; 13(5):605-13.
2. Schmidtchen A *et al.* Acta Derm Venereol 2001; 81(6): 406-9.
3. World Union of Wound Healing Societies. Principles of best practice: Wound exudate and the role of dressings. A consensus document, 2007.
4. Rennie MY *et al.* Understanding Real-Time Fluorescence Signals from Bacteria and Wound Tissues Observed with the MolecuLight *i:X*. Diagnostics (2019).
5. Raizman R *et al.* Use of a bacterial fluorescence imaging device: wound measurement, bacterial detection and targeted debridement. J Wound Care (2019).
6. DaCosta RS *et al.* Point-of-care autofluorescence imaging for real-time sampling and treatment guidance of bioburden in chronic wounds: first-in-human results. PLoS One (2015).
7. Cole W & Coe S. The Use of an Advanced Fluorescence Imaging System to Target Wound Debridement, Decrease Bioburden, Improve Healing Rates, and Provide Positive Revenues in an Outpatient Wound Care Setting. Presented at SAWC Fall 2019 (Las Vegas, NV, USA).

©2019 MolecuLight[®] Inc. All Rights Reserved. PN 1306 Rev 1.1

The MolecuLight[®] *i:X* Imaging Device is approved by Health Canada for sale in Canada and has CE marking for sale in the European Union. The MolecuLight[®] *i:X* Imaging Device has received FDA clearance.

MolecuLight[®] is a Registered Trademark in Canada, the US, and the EU.



Cyan/white color in this image indicates bacteria ($>10^4$ CFU/g)^{4,5}.